VETSCAN®

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1. Webster D, Bignell AHC, Atwood EC. An assessment on the suitability of bromocresol green for the determination of serum albumin. Clin Chim Acta 1974; 53: 101-108.
2. Bergmeyer HU, et al. Provisional recommendations on IFCC methods of catalytic concentrations of enzymes, Part 2. IFCC method for aspartate aminotransferase. Clin Chem 1977; 23: 887-99.
3. Bergmeyer HU, et al. Provisional recommendations on IFCC methods for the measurement of catalytic concentrations of enzymes. Part 2. Revised IFCC method for aspartate aminotransferase. Clin Chem 1978; 24: 720-1.
4. Sampson EJ, et al. A coupled-enzyme equilibrium method for measuring urea in serum: optimization and evaluation of the AACC study group on urea candidate reference method. Clin Chem 1980; 26: 816-826.
5. Cali JP, et al. A reference method for the determination of total calcium in serum. In: Cooper GR, ed., Selected Methods of Clinical Chemistry, Washington, DC: American Association for Clinical Chemistry; 1977; Vol 8: 3-8.
6. Kessler G, Wolfman M. An Automated procedure for the simultaneous determination of calcium and phosphorus. Clin Chem 1964; 10: 686-703.
7. Michaylova V, Ilkova P. Photometric determination of micro amounts of calcium with arsenazo III. Anal Chim Acta 1971; 53: 194-8.
8. Scarpa A, et al. Metallochromic indicators of ionized calcium. Ann NY Acad Sci 1978; 307: 86-112
9. Kuby SA, Noda L, Lardy HA, Adenosinetriphosphate-Creatine Transphosphorylase. J. Biol Chem 1954; 209: 191-201.
10. Tanzer MI, Gilvarg C, Creatine and creatine kinase measurement. J Biol Chem 1959; 01-3204.
11. Nuttall FQ, Wedin DS, Simple rapid colorimetric method for determination of creatine kinase activity. J Lab Clin Med 1966; 24: 332.
12. Oliver IT, A spectrophotometric method for the determination of creatine phosphokinase and myokinase. Biochem 1955; J 61: 116-122.
13. Rosalki SB, An improved procedure or serum creatine phosphokinase determination, J Lab Clin Med 1967; 69: 696-705.
14. Szasz G, Gruber W, Bernt E, Creatine kinase in serum: I. Determination of optimum reaction conditions. Clin Chem 1976; 22: 650-656.
15. Expert panel on enzymes, committee of standards (IFCC). Approval recommendations of IFCC methods for the measurement of catalytic concentrations of enzymes, Part 1. General considerations. Clin Chim Acta, IFCC sections 1979; 98: 163-174.
16. Committee on enzymes of the Scandinavian society for clinical chemistry and clinical physiology. Recommended method for the determination of creatine kinase in blood. Scand J. Clin Lab Invest 1976; 36: 711-723.
17. Knoll VE, et al. Spezifische kreatininbestimmung im serum. Z Klin Chemi Clin Biochem. 1970; 8:582-587.
18. Haeckel R, et al. Simplified determinations of the "true" creatinine concentration in serum and urine. J Clin Chem Clin Biochem 1980; 18: 385-394.
19. Moss GA, Bondar RJL, Buzzelli DM, Kinetic enzymatic method for determining serum creatinine. Clin Chem 1975; 21: 1422-1426.
20. Jaynes PK, Feld RD, Johnson GF. An enzymatic, reaction-rate assay for serum creatinine with a centrifugal analyzer. Clin Chem 1982; 28: 114-117.

21. Fossati P, Prencipe L, Berti G. Enzymatic creatinine assay: A new colorimetric method based on hydrogen peroxide measurement. *Clin Chem* 1983; 29: 1494-1496.
22. Whelton A, et al. Nitrogen Metabolites and Renal Function. In: Tietz Textbook of Clinical Chemistry, 3rd Ed. Burtis CA, Ashwood ER, Eds. Philadelphia: W.B. Saunders Company. 1999;1513-1575.
23. Ball EG, Revel JP, Cooper O. The quantitative measurement of γ -glutamyl transpeptidase activity. *J Biol Chem* 1956; 221: 895-908.
24. Goldbarg JA, et al. The colorimetric determination of γ -glutamyl transpeptidase with a synthetic substrate. *Arch Biochem Biophys*. 1960; 91: 61-70.
25. Orłowski M and Meister A. 1. γ -glutamyl-p-nitroanilide: A new convenient substrate for determination and study of k- and d- γ -glutamyltranspeptidase activities. *Biochem Biophys Acta* 1963; 73: 679-681.
26. Persijn JP, van der Slik W. A new method for the determination of γ -glutamyl-transferase in serum. *J Clin Chem Clin Biochem* 1976; 14: 421-427.
27. Shaw LM et al. IFCC methods for the measurement of catalytic concentration of enzymes. Part 4 IFCC method for γ -glutamyl-transferase. *J Clin Chem Clin Biochem* 1983; 21: 633-646.
28. Somogyi M. A reagent for the copper-idiometric determination of very small amounts of sugar. *J Biol Chem* 1937; 117: 771-776.
29. Nelson N. A photometric adaptation of the Somogyi method for the determination of glucose. *J Biol* 1944; 153: 375-380.
30. Kaplan LA. Glucose. In: Kaplan LA, Pesce AJ, eds. *Clinical Chemistry: Theory, analysis, and correlation*, 2nd ed. St. Louis: The C.V. Mosby Company; 1989; 850-856.
31. Overfield CV, Savory J, Heintges MG. Glycosis: A re-evaluation of the effect on blood glucose. *Clin Chim Acta*. 1972; 39: 35-40.
32. Malloy HT, Evelyn KA. The determination of bilirubin with the photoelectric colorimeter. *J Biol Chem* 1937; 119: 481-90.
33. Meites S. Bilirubin, direct reacting and total, modified Mally-Evelyn method. In: Faulkner WR, Meites S, eds. *Selected Methods of Clinical Chemistry*, Washington, DC: American Association for Clinical Chemistry; 1982; 9: 119-24.
34. Murao S and Tanaka N. A new enzyme "bilirubin oxidase" produced by *Myrothecium verrucaria* MT-1. *Agric Biol Chem* 1981; 45: 2383-4.
35. Osaki S and Anderson S. Enzymatic determination of bilirubin. *Clin Chem* 1982; 30: 971. (Abstract)
36. Perry B, et al. Measurement of total bilirubin by use of bilirubin oxidase. *Clin Chem* 1986; 32: 329-32.
37. Koller A, Kaplan LA. Total serum protein. In: Kaplan LA, Pesce AJ, eds., *Clinical Chemistry: Theory, analysis, and correlation*, 2nd ed. St. Louis: The C.V. Mosby Company; 1989: 1057-60.
38. Reigler E. Eine kolorimetrische bestimmungsmethods des eiweisses. *Z Anal Chem* 1914; 53: 242-5.
39. Weichselbaum TE. An accurate and rapid method for determination of proteins in small amounts of blood, serum, and plasma. *Am J Clin Path* 1946;16: 40-9.
40. Doumas BT, et al. A candidate reference method for determination of total protein in serum. I. Development and validation. *Clin Chem* 1981; 27: 1642-50.
41. National Committee for Clinical Laboratory Standards (NCCLS). Procedures for handling and processing of blood specimens; tentative standard. NCCLS document H18-A2. Wayne, PA: NCCLS, 1999.
42. Overfield CV, Savory J, Heintges MG. Glycosis: A re-evaluation of the effect on blood glucose. *Clin Chim Acta* 1972; 39: 35-40.
43. Rehak NN, Chiang BT. Storage of whole blood: Effect of temperature on the measured concentration of analytes in serum. *Clin Chem* 1988; 34: 2111-14.
44. Sherwin JE, Obernolte R. Bilirubin. In: Kaplan LA, Pesce AJ, eds. *Clinical Chemistry: Theory, Analysis and Correlation*, 2nd ed. St. Louis: The C.V. Mosby Company; 1989; 1009-1015.
45. Henry RJ, Canon DC and Winkelman. *Clinical chemistry principles and technics*, 2nd ed. New York: Harper and Row; 1974; 417-21; 127-8.

46. Witte DL, Brown LF, Williams RL. Effects of bilirubin on detection of hydrogen peroxide by use of peroxidase. Clin Chem 1978; 24: 1778-82.
47. Melnik J, Potter JL. Variance in capillary and venous glucose levels during glucose tolerance test. Am J Med Tech 1982; 48: 543-5.
48. Moss DW, Henderson AR. Enzymes. In: Burtis CA, Ashwood ER, eds. Tietz Textbook of Clinical Chemistry, 2nd edition. Philadelphia: WB Saunders Company. 1994; 804.
49. Van Pelt J. Enzymatic determination of sodium, potassium and chloride in serum compared with determination by flame photometry, coulometry and ion selective electrodes. Clin Chem 1994; 40: 846-847.
50. Hubl W, et al. Enzymatic determination of sodium, potassium and chloride in abnormal (hemolyzed, icteric, lipemic, paraproteinemic, or uremic) serum samples compared with indirect determination with ion selective electrodes. Clin Chem 1994; 40: 1528-1531.
51. Helgersson RC, et al. Host-guest Complexation. 50. Potassium and sodium ion-selective chromogenic ionophores. J Amer Chem Soc 1989; 111: 6339-6350.
52. Kumar A, et al. Chromogenic ionophore-based methods for spectrophotometric assay of sodium and potassium in serum and plasma. Clin Chem 1988; 34: 1709-1712.
53. Berry MN, et al. Enzymatic determination of sodium in serum. Clin Chem 1988; 34: 2295-2298.
54. Malloy HT and Evelyn KA. The determination of bilirubin with the photoelectric colorimeter. J Biol Chem 1937; 119: 481-490.
55. Korzun WJ, Miller WG. Carbon Dioxide. In: Clinical chemistry theory, analysis and correlation, 2nd ed. Kaplan LA, Pesce AJ, eds. St. Louis: The CV Mosby Company. 1989: 869-872.
56. Koller A, Kaplan LA. Total serum protein. In: Clinical Chemistry: Theory, Analysis, and Correlation, 2nd ed. Kaplan LA, Pesce AJ, eds. St. Louis: The C.V. Mosby Company. 1989: 1057-1060.
59. Clinical and Laboratory Standards Institute (formerly, National Committee for Clinical Laboratory Standards, NCCLS). Interference testing in clinical chemistry; proposed guideline. NCCLS Document EP7-A. Wayne, PA: NCCLS, 2002.

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© 2023, Abaxis, Inc., Union City, CA 94587