VETSCAN®

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1. Wróblewski F, LaDue JS. Serum glutamic-pyruvic transaminase in cardiac and hepatic disease. Proc Soc Exp Biol Med 1956; 91: 569-71.
2. Bergmeyer HU, Horder M. IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC Method for alanine aminotransferase. J. Clin Chem Clin Biochem 1980; 18: 521-34.
3. Webster D, et al. An assessment on the suitability of bromocresol green for the determination of serum albumin. Clin Chim Acta 1974; 53:101-8.
4. Bowers GN, et al. IFCC methods for the measurement of catalytic concentration of enzymes. Part 1. General considerations concerning the determination of the catalytic concentration of an enzyme in the blood serum or plasma of man. Clin Chim Acta 1979; 98:163F-74F.
5. Ohmori Y. Über die Phosphomonoesterase. Enzymologia 1937; 4:217-31.
6. Fujita H. Über die Mikrobestimmung der Blutphosphatase. J Biochem, Japan. 1937; 30: 69-87.
7. Petittler C, et al. Mechanism of action of Mg²⁺ and Zn²⁺ on rat placental alkaline phosphatase. I. Studies on the soluble Zn²⁺ and Mg²⁺ alkaline phosphatase. Can J Biochem 1975; 53: 1089-1100.
8. Tietz NW, et al. A reference method for measurement of alkaline phosphatase activity in human serum. Clin Chem 1983; 29: 751-61.
9. Bergmeyer HU, et al. Provisional recommendations on IFCC methods of catalytic concentrations of enzymes, Part 2. IFCC method for aspartate aminotransferase. Clin Chem 1977; 23: 887-99.
10. Bergmeyer HU, et al. Provisional recommendations on IFCC methods for the measurement of catalytic concentrations of enzymes. Part 2. Revised IFCC method for aspartate aminotransferase. Clin Chem 1978; 24: 720-1.
11. Fales FW. Urea in serum, direct diacetyl monoxime method. In: Selected Methods of Clinical Chemistry, Vol 9. Faulkner WR, Meites S, eds. Washington, DC.: American Association for Clinical Chemistry. 1982: 365-373.

12. Van Slyke, et al. A permanent preparation of urease, and its use in the determination of urea. *J Biol Chem* 1914; 19: 11-228.
13. Fawcett JK, et al. A rapid and precise method for the determination of urea. *J Clin Pathol* 1960; 13: 156-159.
14. Chaney, et al. Urea and ammonia determinations. *Clin Chem* 1962; 8: 130-132.
15. Talke H, et al. Enzymatische harnstoffbestimmung in blut und serum im optischen Test nach Warburg. *Klin Wochensh*, 1965; 43: 174-175.
16. Hallett, et al. Reduced nicotinamide adenine dinucleotide-coupled reaction for emergency blood urea estimation. *Clin Chim Acta* 1971; 35: 33-37.
17. Patton, et al. Spectrophotometric and kinetics investigation of the Berthelot reaction for the determination of ammonia. *Anal Chem* 1977; 49: 464-469.
18. Sampson EJ, et al. A coupled-enzyme equilibrium method for the measuring urea in serum: optimization and evaluation of the AACC study group on Urea Candidate reference method. *Clin Chem* 1980; 26: 816-826.
19. Cali JP, Bowers GN, Young DS, et al. A reference method for the determination of total calcium in serum. In: GR Cooper, ed., *Selected methods of Clinical Chemistry*. Washington, DC: American Association for Clinical Chemistry 1977; 8: 3-8.
20. Kessler G, M Wolfman. An Automated procedure for the simultaneous determination of calcium and phosphorus. *Clin Chem* 1964; 10: 686-703.
21. Michaylova V, P Ilkova. Photometric determination of micro amounts of calcium with arsenazo III. *Anal Chim Acta* 1971; 53: 194-8.
22. Scarpa A, et al. Metallochromic indicators of ionized calcium. *Ann NY Acad Sci* 1978; 307: 86-112.
23. Ono T, et al. A new enzymatic assay of chloride in serum. *Clin Chem* 1988; 34: 552-553.
24. Knoll VE, et al. Spezifische kreatininbestimmung Im Serum. *Z Klin Chemi Clin Biochem* 1970; 8: 582-587.
25. Haeckel R, et al. Simplified determinations of the "true" creatinine concentration in serum and urine. *J Clin Chem Clin Biochem* 1980; 18: 385-394.
26. Moss GA, et al. Kinetic enzymatic method for determining serum creatinine. *Clin Chem* 1975; 21: 1422-1426.
27. Jaynes PK, et al. An enzymatic, reaction-rate assay for serum creatinine with a centrifugal analyzer. *Clin Chem* 1982; 28: 114-117.
28. Fossati P, et al. Enzymatic creatinine assay: a new colorimetric method based on hydrogen peroxide measurement. *Clin Chem* 1983; 29: 1494-1496.
29. Whelton A, et al. Nitrogen metabolites and renal function. In: CA Burtis and ER Ashwood, Eds., *Tietz Textbook of Clinical Chemistry*, 3rd Ed. Philadelphia: W.B. Saunders Company. 1999; 1513-1575.
30. Folin O, and Wu H. A system of blood analysis. *J Biol Chem* 1919; 38: 81-110.
31. Somogyi M. A reagent for the copper-idiometric determination of very small amounts of sugar. *J Biol Chem* 1937; 117: 771-776.
32. Nelson N. A photometric adaption of the Somogyi method for the determination of glucose. *J Biol* 1944; 153: 375-380.
33. Kaplan LA. Glucose. In: LA Kaplan and AJ Pesce, eds., *Clinical Chemistry: Theory, Analysis, and Correlation*, 2nd ed St. Louis: The C.V. Mosby Company. 1989; 850-856.
34. Berry MN, et al. Enzymatic determination of potassium in serum. *Clin Chem* 1989; 35: 817-20.
35. Van Pelt J. Enzymatic determination of sodium, potassium and chloride in serum compared with determination by flame photometry, coulometry and ion selective electrodes. *Clin Chem* 1994; 40: 846-7.
36. Hubl W, et al. Enzymatic determination of sodium, potassium and chloride in abnormal (hemolyzed, icteric, lipemic, paraproteinemic, or uremic) serum samples compared with indirect determination with ion selective electrodes. *Clin Chem* 1994; 40: 1528-31.

37. Helgerson RC, et al. Host-guest Complexation. 50. Potassium and sodium ion-selective chromogenic ionophores. *J Amer Chem Soc* 1989; 111: 6339-50.
38. Kumar A, et al. Chromogenic ionophore-based methods for spectrophotometric assay of sodium and potassium in serum and plasma. *Clin Chem* 1988; 34: 1709-12.
39. Berry MN, et al. Enzymatic determination of sodium in serum. *Clin Chem* 1988; 34: 2295-8.
40. Malloy HT, Evelyn KA. The determination of bilirubin with the photoelectric colorimeter. *J Biol Chem* 1937; 119: 481-90.
41. Meites S. Bilirubin, direct reacting and total, modified Mally-Evelyn method. In: WR Faulkner and S Meites, eds., *Selected Methods of Clinical Chemistry*, Washington, DC: American Association for Clinical Chemistry. 1982; 9: 119-24.
42. Murao S, Tanaka N. A new enzyme "bilirubin oxidase" produced by *Myrothecium verrucaria* MT-1. *Agric Biol Chem* 1981; 45: 2383-4.
43. Osaki S, Anderson S. Enzymatic determination of bilirubin. *Clin Chem* 1982; 30: 971 (Abstract).
44. Perry B, et al. Measurement of total bilirubin by use of bilirubin oxidase. *Clin Chem* 1986; 32: 329-32.
45. Skeggs LT Jr. An automatic method for the determination of carbon dioxide in blood plasma. *Am J Clin Pathol* 1960; 33: 181-185.
46. Korzun WJ, Miller WG. Carbon Dioxide. In: *Clinical chemistry theory, analysis and correlation*, 2nd ed. Kaplan LA, Pesce AJ, eds. St. Louis: The CV Mosby Company. 1989: 869-872.
47. Koller A, Kaplan LA. Total serum protein. In: LA Kaplan and AJ Pesce, eds., *Clinical Chemistry: Theory, Analysis, and Correlation*, 2nd ed. St. Louis: The C.V. Mosby Company. 1989: 1057-60.
48. Reigler E. Eine kolorimetrische Bestimmungsmethods des Eiweisses. *Z Anal Chem* 1914; 53: 242-5.
49. Weichselbaum TE. An accurate and rapid method for determination of proteins in small amounts of blood serum and plasma. *Am J Clin Path* 1946; 16: 40-9.
50. Doumas BT, et al. A candidate reference method for determination of total protein in serum. I. Development and validation. *Clin Chem* 1981; 27: 1642-50.
51. CLSI. *Procedures for Handling and Processing of Blood Specimens*; tentative standard. CLSI document H18-A2. Wayne, PA: CLSI, 1999.
52. Overfield CV, Savory J, and Heintges MG. Glycosis: a re-evaluation of the effect on blood glucose. *Clin Chim Acta* 1972; 39: 35-40.
53. Rehak NN, Chiang BT. Storage of whole blood: effect of temperature on the measured concentration of analytes in serum. *Clin Chem* 1988; 34: 2111-14.
54. Sherwin JE, Obernolte R. Bilirubin. In: LA Kaplan and AJ Pesce, eds. *Clinical Chemistry: Theory, Analysis and Correlation*, 2nd ed. St. Louis: The C.V. Mosby Company. 1989; 1009-1015.
55. Henry RJ, Canon DC and Winkelman. *Clinical Chemistry Principles and Technics*, 2nd ed. New York: Harper and Row. 1974; 417-421: 127-8.
56. Scott MG. Electrolytes and Blood Gases. In: *Tietz Textbook of Clinical Chemistry*. 3rd ed. Burtis CA, Ashwood ER, eds. Philadelphia: WB Saunders Company. 1999: 1065-1066.
57. Scott MG, Electrolytes and Blood Gases. In: Burtis CA, Ashwood ER, eds. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: WB Saunders Company. 1999: 617-721.
58. Witte DL, Brown LF, and RL Williams. Effects of bilirubin on detection of hydrogen peroxide by use of peroxidase. *Clin Chem* 1978; 24: 1778-82.
59. Melnik J, Potter JL. Variance in capillary and venous glucose levels during glucose tolerance test. *Am J Med Tech* 1982; 48: 543-5.
60. Henry RJ, Canon DC and Winkelman. *Clinical chemistry principles and technics*, 2nd ed. New York: Harper and Row. 1974; 417-21: 127-8.

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© 2023, Abaxis, Inc., Union City, CA 94587