VetScan®

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1. Wróbleski F and LaDue. Serum glutamic-pyruvic transaminase in cardiac and hepatic disease. Proc Soc Exp Biol Med. 1956;91:569-71.
2. Bergmeyer HU and Horder M. IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase. J. Clin Chem Clin Biochem 1980;18:521-34.
3. Webster D, et al. An assessment on the suitability of bromocresol green for the determination of serum albumin. Clin Chim Acta 1974;53:101-8.
4. Bowers GN, et al. IFCC methods for the measurement of catalytic concentration of enzymes. Part 1. General considerations concerning the determination of the catalytic concentration of an enzyme in the blood serum or plasma of man. Clin Chim Acta 1979;98:163F-74F.
5. Ohmori Y. Über die Phosphomomesterase. Enzymologia 1937;4:217-31.
6. Fujita H. Über die Mikrobestimmung der Blutphosphatase. J Biochem, Japan. 1937;30:69-87.
7. Petitclerc C, et al. Mechanism of action of Mg^{2+} and Zn^{2+} on rat placental alkaline phosphatase. I. Studies on the soluble Zn^{2+} and Mg^{2+} alkaline phosphatase. Can J Biochem 1975;53:1089-1100.
8. Tietz NW, et al. A reference method for measurement of alkaline phosphatase activity in human serum. Clin Chem 1983;29:751-61.
9. Somogyi M. Modifications of two methods for the assay of amylase. Clin Chem 1960;6:23-35.
10. Gillard BK, Markman HC and Feig SA. Direct spectrophotometric determination of α -amylase activity in saliva with p-nitrophenyl α -maltoside as substrate. Clin Chem 1977;23:2279-82.
11. Wallenfels K, et al. The enzymic synthesis, by transglucosylation of a homologous series of glycosidically substituted malto-ologosaccharides, and their use as amylase substrates. Carbohydrate Res 1978;61:359-68.
12. Cali JP, Bowers GN, Young DS, et al. A reference method for the determination of total calcium in serum. In: GR Cooper, ed., Selected methods of Clinical Chemistry. Washington, DC: American Association for Clinical Chemistry; 1977;8:3-8.
13. Kessler G, M Wolfman. An Automated procedure for the simultaneous determination of calcium and phosphorus. Clin Chem 1964;10: 686-703.
14. Michaylova V, P Ilkova. Photometric determination of micro amounts of calcium with arsenazo III. Anal Chim Acta 1971;53: 194-8.
15. Scarpa A, et al. Metallochromic indicators of ionized calcium. Ann NY Acad Sci 1978;307:86-112.
16. Knoll VE, et al. Spezifische Kreatininbetimmung Im Serum. Z Klin Chemi Clin Biochem. 1970;8:582-587.
17. Haeckel R, et al. Simplified Determinations of the "True" Creatinine Concentration In Serum And Urine. J Cklin Chem Clin Biochem. 1980;18:385-394.
18. Moss GA, et al. Kinetic Enzymatic Method For Determining Serum Creatinine. Clin Chem 1975;21:1422-1426.
19. Jaynes PK, et al. An Enzymatic, Reaction-Rate Assay For Serum Creatinine With a Centrifugal Analyzer. Clin Chem 1982;28:114-117.
20. Fossati P, et al. Enzymatic Creatinine Assay: A New Colorimetric Method Based on Hydrogen Peroxide Measurement. Clin Chem 1983;29:1494-1496.
21. Whelton A, et al. Nitrogen Metabolites and Renal Function. In:CA Burtis and ER Ashwood, Eds., Tietz Textbook of Clinical Chemistry, 3rd Ed. Philadelphia: W.B. Saunders Company. 1999;1513-1575.

22. Folin O, and Wu H. A System of blood analysis. J Biol Chem 1919;38: 81-110.
23. Somogyi M. A reagent for the copper-idiometric determination of very small amounts of sugar. J Biol Chem 1937;117: 771-776.
24. Nelson N. A photometric adaption of the Somogyi method for the determination of glucose. J Biol 1944; 153: 375-380.
25. Kaplan LA. Glucose. In: LA Kaplan and AJ Pesce, eds., Clinical Chemistry: Theory, Analysis, and Correlation, 2nd ed St. Louis: The C.V. Mosby Company; 1989;850-856.
26. Schulz DW, et al. An Enzymic Method for the Measurement of Inorganic Phosphate Determination Anal Biochem 1967;19:300-14.
27. Tedokon, M Suzuki, et al. Enzymatic Assay of Inorganic Phosphate with Use of Sucrose Phosphorylase and Phosphoglucomutase. Clin Chem 1992;38:512-5.
28. Berry MN, et al. Enzymatic determination of potassium in serum. Clin Chem 1989;35:817-20.
29. Van Pelt J. Enzymatic determination of sodium, potassium and chloride in serum compared with determination by flame photometry, coulometry and ion selective electrodes. Clin Chem 1994;40:846-7.
30. Hubl W, et al. Enzymatic determination of sodium, potassium and chloride in abnormal (hemolyzed, icteric, lipemic, paraproteinemic, or uremic) serum samples compared with indirect determination with ion selective electrodes. Clin Chem 1994;40:1528-31.
31. Helgersson RC, et al. Host-guest Complexation. 50. Potassium and sodium ion-selective chromogenic ionophores. J Amer Chem Soc 1989;111:6339-50.
32. Kumar A, et al. Chromogenic ionophore-based methods for spectrophotometric assay of sodium and potassium in serum and plasma. Clin Chem 1988;34:1709-12.
33. Berry MN, et al. Enzymatic determination of sodium in serum. Clin Chem 1988;34:2295-8.
34. Malloy HT, and Evelyn KA. The determination of bilirubin with the photoelectric colorimeter. J Biol Chem 1937; 119:481-90.
35. Meites S. Bilirubin, direct reacting and total, modified Mally-Evelyn method. In: WR Faulkner and S Meites, eds., Selected Methods of Clinical Chemistry, Washington, DC: American Association for Clinical Chemistry; 1982;9:119-24.
36. Murao S and Tanaka N. A new enzyme "bilirubin oxidase" produced by Myrothecium verrucaria MT-1. Agric Biol Chem 1981;45:2383-4.
37. Osaki S, and S Anderson. Enzymatic determination of bilirubin. Clin Chem 1982;30:971. (Abstract)
38. Perry B, et al. Measurement of total bilirubin by use of bilirubin oxidase. Clin Chem 1986;32:329-32.
39. Koller A and Kaplan LA. Total serum protein. In: LA Kaplan and AJ Pesce, eds., Clinical Chemistry: Theory, Analysis, and Correlation, 2nd ed. St. Louis: The C.V. Mosby Company; 1989:1057-60.
40. Reigler E. Eine kolorimetrische Bestimmungsmethods des Eiweisses. Z Anal Chem 1914;53:242-5.
41. Weichselbaum TE. An accurate and rapid method for determination of proteins in small amounts of blood serum and plasma. Am J Clin Path 1946;16:40-9.
42. Doumas BT, et al. A candidate reference method for determination of total protein in serum. I. Development and validation. Clin Chem 1981;27:1642-50.
43. Sampson, EJ MA Baird, CA Burtis, EM Smith, DL Witte, and DD Bayse. A coupled-enzyme equilibrium method for measuring urea in serum: optimization and evaluation of the AACC study group on urea candidate reference method. Clin Chem 1980;26: 816-826.
44. National Committee for Clinical Laboratory Standards (NCCLS). Procedures for Handling and Processing of Blood Specimens; tentative standard. NCCLS document H18-A2. Wayne, PA: NCCLS, 1999.

45. Overfield CV, Savory J, and Heintges MG. Glycosis: a re-evaluation of the effect on blood glucose. *Clin Chim Acta* 1972;39:35-40.
46. Rehak NN and Chiang BT. Storage of whole blood: effect of temperature on the measured concentration of analytes in serum. *Clin Chem* 1988;34:2111-14.
47. Sherwin JE and Oberholte R. Bilirubin. In: LA Kaplan and AJ Pesce, eds. *Clinical Chemistry: Theory, Analysis and Correlation*, 2nd ed. St. Louis: The C.V. Mosby Company; 1989;1009-1015.
48. Henry RJ, Canon DC and Winkelman. *Clinical Chemistry Principles and Technics*, 2nd ed. New York: Harper and Row; 1974;417-21; 127-8.
49. Scott MG, Electrolytes and Blood Gases. In: Burtis CA, Ashwood ER, eds. *Tietz Textbook of Clinical Chemistry*. 3rd ed. Philadelphia: WB Saunders Company, 1999;617-721.
50. Witte DL, Brown LF, and RL Williams. Effects of bilirubin on detection of hydrogen peroxide by use of peroxidase. *Clin Chem* 1978;24:1778-82.
51. Melnik J and Potter JL. Variance in capillary and venous glucose levels during glucose tolerance test. *Am J Med Tech* 1982;48:543-5.

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© 2023, Abaxis, Inc., Юніон-Сіті, Каліфорнія, 94587